

iSWAB – Beckman Coulter Biomek Compatibility Guide

Biomek NX/FX Span 8 Head

Tips: Span-8 P20, P50 or P250 Tips

<https://www.beckmancoulter.com/wsrportal/wsr/research-and-discovery/products-and-services/research-automation/tips/index.htm>

Rack holder: 24-Position Tube Rack (Part # 373661) with 11 mm white tube insert (Part # 373696)

<https://www.beckmancoulter.com/wsrportal/page/itemDetails?itemNumber=373661#2/10//0/25/1/0/asc/2/373661///0/1//0/>



Chemistry: Genfind v2 and/or DNAdvance

